



## **Pepha®-Tight**

**Pepha®-Tight** is a highly purified biotechnologically produced extract from the microalgae *Nannochloropsis oculata* combined with a well balanced fraction of polysaccharides. It functions as a skin tightener which provides immediate tightening to the skin. It has been shown to protect human fibroblasts from oxidative stress, simultaneously increasing the formation of collagen-I thereby promoting a long term tightening effect through the strengthening of the skin's connective tissue.

**Pepha®-Tight** can be used for all skin tightening applications including lifting products, masks, decollete treatments, gels and serums.

**Usage Rate:** 1 - 5%

**Appearance:** Yellow to amber, opalescent to turbid, highly viscous liquid

**pH:** 5.1

**Solubility:** Soluble in water.

**Storage:** Clean location, protected from light, temperature between 35° - 46°F (2°-8°C)

**INCI:** Algae Extract (and) Pullulan

**Caution: Do not directly apply onto your skin.**

This product should be added to a formulation at the recommended usage rate.

### **Disclaimer**

Statements contained herein have not been evaluated by the Food and Drug Administration. These products are not intended to diagnose, treat and cure or prevent disease. Always consult with your professional skin care provider.

### **Based on 1 ounce formulation:**

**1% of 1 oz – 1/8 tsp** to 1 oz water/mixture

**2% of 1 oz – 1/4 tsp** to 1 oz water/mixture

**4% of 1 oz – 1/2 tsp** to 1 oz water/mixture

**5% of 1 oz – 2/3 tsp** to 1 oz water/mixture

### **For a 3% formulation:**

To a **4 oz** base for **3%** add **1 1/4 tsp**

To a **1-gallon** base for **3%** add **3/4 cup**

**REORDER CONTACT:** Anita Nelson @ModelSupplies on Twitter & Facebook~!

eBay: <http://Model-Supplies.com>

Etsy: <http://eModelSupplies.com>

Blog: <http://ModelSuppliesBlog.com>

Radio: <http://ModelSuppliesRadio.com>

Video: <http://YouTube.com/ModelSupplies>